

#### QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE			
<b>Department:</b> Quality Control	SOP No.:		
Title: Handling and Usage of Reagents and Chemicals	<b>Effective Date:</b>		
Supersedes: Nil	Review Date:		
Issue Date:	Page No.:		

#### 1.0 OBJECTIVE:

To lay down a procedure for handling & usage of reagents & chemicals.

#### **2.0 SCOPE**:

This SOP is applicable to handling & usage of reagents & chemicals in Quality Control Department.

#### 3.0 **RESPONSIBILITY:**

Officer, Executive – Quality Control Department Head – Quality Control Department

#### 4.0 **DEFINITION**(S):

NA

#### **5.0 PROCEDURE:**

- 5.1 After receipt of Reagents and Chemicals from suppliers, check the Grade/Class of material like Reagents / Solvents items against the details given in delivery challan and verify whether same Reagent /chemical was indented.
- 5.2 Check the integrity of pack before placing on rack.
- 5.3 Arrange the chemicals in dedicated place in the quality control laboratory and affix the received labels as per Annexure-I.
- 5.4 Receipt of reagents, solvent and acids and issuance records should be recorded as per Annexure II.
- Assign shelf life 5 years, whenever expiry is not mentioned. Maintain the date of opening on Chemical /Reagent, Solvent bottle at the time of opening for use and also mention use before date on the label as per Annexure-I. All opened reagent and chemicals shall be used within 2 years from the date of opening or expiry date declared by the manufacturer, which ever is earlier.
- 5.6 Reagents shall be visually examined during their usage.
- 5.7 If any discoloration, lump formation or precipitation (in case of liquids) or any unusual appearance is observed, it shall not be used and shall be discarded.



QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE		
<b>Department:</b> Quality Control	SOP No.:	
Title: Handling and Usage of Reagents and Chemicals	<b>Effective Date:</b>	
Supersedes: Nil	<b>Review Date:</b>	
Issue Date:	Page No.:	

- 5.8 All reagents and Chemicals shall be stored at the room temperature, unless any other storage condition is specified on the label.
- 5.9 Toxic chemicals / reagents shall be identified.
- 5.10 All toxic chemicals / reagents shall be stored separately under lock and key.
- 5.11 List of all toxic chemicals shall be maintained in quality control department.
- 5.12 MSDS (Material Safety Data Sheets) of all toxic chemicals shall be maintained in quality control department.
- 5.13 Maintain the destruction record of all reagents / solvents after the expiry date or rejection due to any discoloration, lump formation or precipitate, as per Annexure –III

#### **6.0 ABBREVIATION** (S):

QCD - Quality Control Department

**SOP - Standard Operating Procedure** 

#### 7.0 REFERENCE (S):

NA

#### 8.0 ANNEXURE (S):

Annexure -I Receipt & opening label

Annexure-II Stock Record of Reagents, Solvents and Acids

Annexure -III Destruction record of Reagents/Chemicals

#### 9.0 REVISION CARD:

S.No.	REVISION No.	REVISION DATE	DETAILS OF REVISION	REASON (S) FOR REVISION		



QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE				
Department: Quality Control SOP No.:				
Title: Handling and Usage of Reagents and Chemicals	Effective Date:			
Supersedes: Nil	Review Date:			
Issue Date:	Page No.:			

# ANNEXURE I Receipt & Opening label

Reagent			
Date of Receipt	1		
Assigned Shelf life / Expiry date	2		
Date of Opening	3		
Use before	4		
Sign. Date	5		



QUALITY CONTROL DEPARTMENT

STANDARD OPERATING PROCEDURE				
Department: Quality Control	SOP No.:			
Title: Handling and Usage of Reagents and Chemicals	Effective Date:			
Supersedes: Nil	Review Date:			
Issue Date:	Page No.:			

### ANNEXURE II Stock Record of Reagents, Solvents and Acids

S.No	Date of Received	Opening Stock	Receiving Qty.	Issue Qty.	Closing Record	Date of Issue	Issue By	Remarks