

PHARMA DEVILS ENGINEERING DEPARTMENT

STANDARD OPERATING PROCEDURE	
Department: Engineering	SOP No.:
Title: Specification of Steam Filter	Effective Date:
Supersedes: Nil	Review Date:
Issue Date:	Page No.:

1.0 **Revision History**

Rev. No.	Details of changes	Reason for change
00	NIL	NEW SOP

2.0 **OBJECTIVE:**

2.1 The Objective of this SOP is to provide the Specifications of Steam Filter.

3.0 **SCOPE:**

3.1 This SOP is applicable for the provide the Specifications of Steam Filters.

4.0 **RESPONSIBILITY:**

4.1 The Maintenance Engineer shall be:

4.1.1 Responsible for ensuring the specification at the time of releasing of new order.

4.1.2 Responsible to ensure the compliance of specifications at the time of receipt.

5.0 **ACCOUNTABILITY:**

Head – Engineering Services

PROCEDURE: 6.0

6.1 Following are the specifications of the filter cartridge being used for steam.

6.2 **Final Filters:**

- 6.2.1 Cartridge Porosity
- 6.2.2 Material of construction
- 6.2.3 Application
- 6.2.4 Manufacturer
- 6.2.5 Supplier

- : 1µ (Micron)
- : Sintered SS 316L
- : Pure Steam and Raw Steam
- : M/S Pall filters.
- : Pall Pharmalab (India)



PHARMA DEVILS

ENGINEERING DEPARTMENT

STANDARD OPERATING PROCEDUREDepartment: EngineeringSOP No.:Title: Specification of Steam FilterEffective Date:Supersedes: NilReview Date:Issue Date:Page No.:

7.0 ANNEXURES:

NIL

8.0 References (S)

NIL

9.0 Glossary

SOP	: Standard Operating procedure
No	: Number