



# PHARMA DEVILS

## ENGINEERING DEPARTMENT

### STANDARD OPERATING PROCEDURE

<b>Department:</b> Engineering	<b>SOP No.:</b>
<b>Title:</b> Specification of Vent Filters	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

#### 1.0 Revision History

Rev. No.	Details of changes	Reason for change
00	NIL	NEW SOP

#### 2.0 OBJECTIVE:

2.1 The Objective of this SOP is to provide the Specifications of Vent Filters.

#### 3.0 SCOPE:

3.1 This SOP is applicable for the provide the Specifications of Vent Filters.

#### 4.0 RESPONSIBILITY:

4.1 The Maintenance Engineer shall be:

4.1.1 Responsible for ensuring the specification at the time of releasing of new order.

4.1.2 Responsible to ensure the compliance of specifications at the time of receipt.

#### 5.0 ACCOUNTABILITY:

Head –Engineering Services

#### 6.0 PROCEDURE:

6.1 Following are the specifications of the filter cartridge being used for vent filters.

##### 6.2 Final Filters:

- 6.2.1 Pore size : 0.2  $\mu$  (Micron)
- 6.2.2 Type : Hydrophobic
- 6.2.3 Material of construction : Class VI plastic
- 6.2.4 Manufacturer : M/S Pall filtration.
- 6.2.5 Supplier : Pall Pharmalab (India)
- 6.2.6 Application : Storage tank of PW and WFI Autoclave



# PHARMA DEVILS

## ENGINEERING DEPARTMENT

### STANDARD OPERATING PROCEDURE

<b>Department:</b> Engineering	<b>SOP No.:</b>
<b>Title:</b> Specification of Vent Filters	<b>Effective Date:</b>
<b>Supersedes:</b> Nil	<b>Review Date:</b>
<b>Issue Date:</b>	<b>Page No.:</b>

#### 7.0 ANNEXURES:

NIL

#### 8.0 References (S)

NIL

#### 9.0 Glossary

SOP : Standard Operating procedure

No : Number