



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

ANNEXURE I

AFTER 1st BATCH MANUFACTURING

S.No.	Equipment No.	Location of Sample	Residue Detection
1.	Stirrer Tank,	1. Stirrer	
		2. Inner Surface of Stirrer Tank	
		3. Lid	
2.	Sugar tank,	1. Bottom valve	
		2. Lid	
		3. Stirrer	
		4. Inner surface	
		5. Top dish	
3.	Mfg. tank,	1. Bottom valve	
		2. Lid	
		3. Stirrer	
		4. Inner surface	
		5. Top dish	
4.	Transfer pump,	1. Inlet	
		2. Out let	
5.	Storage tank,	1. Bottom valve	
		2. Lid	
		3. Stirrer	
		4. Inner surface	
		5. Top dish	
6.	Filter press,	1. Filter plate	
		2. Inlet	
		3. Outlet	
7.	Basket filter	1. Inlet	
		2. Outlet	
		3. Screen	
8.	Homogeniser,	1. Inlet	
		2. Discharge chute	
9.	Filling machine,	1. Inlet line	
		2. Filling nozzle	
		3. Storage vessel	
		4. Hose	
		5. Silicon Tube	

Checked By:

Reviewed By:



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

CLEANING VALIDATION PROTOCOL & REPORT FOR MEFENAMIC ACID IP

AFTER IInd BATCH MANUFACTURING

S.No.	Equipment No.	Location of Sample	Residue Detection
1.	Stirrer Tank	1. Stirrer	
		2. Inner Surface of Stirrer Tank	
		3. Lid	
2.	Sugar tank	1. Bottom valve	
		2. Lid	
		3. Stirrer	
		4. Inner surface	
		5. Top dish	
3.	Mfg. tank	1. Bottom valve	
		2. Lid	
		3. Stirrer	
		4. Inner surface	
		5. Top dish	
4.	Transfer pump	3. Inlet	
		4. Out let	
5.	Storage tank	1. Bottom valve	
		2. Lid	
		3. Stirrer	
		4. Inner surface	
		5. Top dish	
6.	Filter press	4. Filter plate	
		5. Inlet	
		6. Outlet	
7.	Basket filter	4. Inlet	
		5. Outlet	
		6. Screen	
8.	Homogeniser	1. Inlet	
		2. Discharge chute	
9.	Filling machine	1. Inlet line	
		2. Filling nozzle	
		3. Storage vessel	
		4. Hose	
		5. Silicon Tube	

Checked By:

Reviewed By: