



PHARMADEVILS

QUALITY ASSURANCE

CLEANING VALIDATION PROTOCOL

Document Name: Cleaning Validation Protocol For Cefixime-100 DT Tablets

Document No.:

Rev No.: 00

Page 1 of 3

Effective Date:

Annexure-6
Visual Inspection Results



PHARMADEVILS

QUALITY ASSURANCE

CLEANING VALIDATION PROTOCOL

Document Name: Cleaning Validation Protocol For Cefixime-100 DT Tablets

Document No.:

Rev No.: 00

Page 2 of 3

Effective Date:

S.No.	Sampling Location	Batch Number of Previous Product			Acceptance limit
1.	Vibro Sifter				
	Outlet of Feeder Hopper				Visually clean
	Inner Surface of Sifter Bowl				
	Sifter Sieve surface				
	Gasket of Sifter Ring				
2.	Multimill				
	Inner surface				Visually clean
	Discharge Valve				
3.	Roll Compactor				
	Hopper				Visually clean
	Roll of Machine				
	Outlet				
4.	Octagonal Blender				
	Side Wall Corner				Visually clean
	Discharge Port				
	Gasket				
5.	Tray Dryer				
	Upper Tray				Visually clean
	Middle Tray				
	Lower Tray				
6.	Compression Machine-16 Stations				
	Turret				Visually clean
	Discharge Chute				
	Inner Surface of Hopper				
	Punch Tip				
7.	Compression Machine-27 Stations				
				Page No.	Page 2 of 3



PHARMADEVILS

QUALITY ASSURANCE

CLEANING VALIDATION PROTOCOL

Document Name: Cleaning Validation Protocol For Cefixime-100 DT Tablets

Document No.:

Rev No.: 00

Page 3 of 3

Effective Date:

S.No.	Sampling Location	Batch Number of Previous Product			Acceptance limit
	Turret				Visually clean
	Discharge Chute				
	Inner Surface of Hopper				
	Punch Tip				
8.	Auto Coater				
	Inner surface of pan				Visually clean
	Edges of pan				
9.	Strip Packing Machine				
	Hopper				Visually clean
	Tablet Falling Chute				
	Bowl				
10.	Blister Packing Machine				
	Hopper				Visually clean
	Tablet Falling Chute				
	Feeding Box				
11.	Alu-Alu Packing Machine				
	Hopper				Visually clean
	Tablet Falling Chute				
	Feeding Box				

Conclusion:

Checked By:
(Signature/Date)

Verified By:
(Signature/Date)

Approved By:
(Signature/Date)