



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

STANDARD TESTING PROCEDURE FOR PRINTED ALUMINIUM FOIL

SAMPLING AND HANDLING

PRESTORAGE REQUIREMENT:

Store in well - closed containers.

SAMPLING:

Carry out the sampling as per SOP No.

QUANTITY OF COMPOSITE SAMPLE FOR ANALYSIS:

30 gm

QUANTITY OF RESERVE SAMPLE:

2 x 30 gm

RETEST PERIOD:

One year

HAZARDS AND PRECAUTIONS, IF ANY:

Use hand gloves and nose mask while sampling.

DISPOSITION OF ANALYTICAL SAMPLE:

Carry out disposition of analytical sample as per SOP No.

TESTING PROCEDURE

DESCRIPTION:

Examine the composite samples and note the observation.

SOLUBILITY:

Insoluble in solutions of alkalis.

IDENTIFICATION:

Determine the melting point using suitable melting point apparatus

RELATED SUBSTANCES (BY TLC):

Plate Type : Silica gel GF 254

Mobile Phase : Mixture of 60 volume of chloroform: 10 volume of ether and 10 volume of glacial acetic acid

Application Volume : 10µl