



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

ANALYTICAL METHOD VALIDATION REPROT FOR LIVOFLOXACIN TABLETS (FOR ASSAY)

S. No	Analytical parameter	Result	Acceptance criteria
1.0	Specificity a) Blank interference: b) Placebo interference:		There shall be no interference in Sample response due to blank and placebo.
2.0	Precision: 2.1 System precision: For		
	Area Replicate	Area Response	RSD NMT 2.0 %.
	1		
	2		
	3		
	4		
	5		
	6		
	Average		
	% RSD		
	System Suitability		
	Tailing factor		Not more than 1.8
	% Standard RSD		RSD NMT 2.0 %.
	2.2 Method precision:		
	Results Replicate	Results (%)	Six results RSD NMT 2.0 %.
	1		
	2		
	3		
	4		
	5		
	6		
	Average		
	% RSD		
	System Suitability		
	Tailing factor		Not more than 1.8
	% Standard RSD		RSD NMT 2.0 %.
	2.3 Intermediate precision: For		



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

ANALYTICAL METHOD VALIDATION REPROT FOR LIVOFLOXACIN TABLETS (FOR ASSAY)

S. No	Analytical parameter	Result	Acceptance criteria	
	Results Replicate	Results (%)	Six results RSD NMT 2.0 %.	
	1			
	2			
	3			
	4			
	5			
	6			
	Average			
	% RSD			
	System Suitability			
	Tailing factor		Not more than 1.8	
	% Standard RSD		RSD NMT 2.0 %.	
	Analyst 1-			
	Batch No.	Components	Assay %	Result RSD NMT 2.0 %.
Analyst 2-				
Batch No.	Components	Assay %		
Result % RSD				
3.0	Linearity			
	Level	Area Response	RSD NMT 2.0 %.	
	80 %			
	90 %			
	100 %			
	110 %			
	120%			
	Average			
	% RSD			
	System Suitability			
	Tailing factor		Not more than 1.8	
	% Standard RSD		RSD NMT 2.0 %.	



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

ANALYTICAL METHOD VALIDATION REPROT FOR LIVOFLOXACIN TABLETS (FOR ASSAY)

S. No	Analytical parameter	Result	Acceptance criteria		
	Coefficient of determination (r2) .		(r2) NMT- 0.999		
4.0	Range				
	Lower Level - For 80 % level				
	Area Replicate	Area Response			
	1				
	2				
	3				
	4				
	5				
	6				
	Average				
	Tailing factor		Not more than 1.8		
	% Standard RSD		RSD NMT 2.0 %.		
	Lower Level - For 120 % level				
	Area Replicate	Area Response			
	1				
	2				
	3				
	4				
	5				
	6				
	Average				
	Tailing factor				
	% Standard RSD				
	Acceptance Criteria				
% RSD of upper concentration		RSD NMT 2.0 %.			
% RSD of lower concentration		RSD NMT 2.0 %.			
Correlation Coefficient		(r2) should be greater than 0.99			
5.0	Accuracy and Recovery				
	S. No	Sample	Recovery	Spike	Average



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

ANALYTICAL METHOD VALIDATION REPROT FOR LIVOFLOXACIN TABLETS (FOR ASSAY)

S. No	Analytical parameter	Result			Acceptance criteria	
		in mg	in mg	results in %	Each individual sample recovery should lie within the range of 98% to 102%.	
1.	80 % sample					
2.	100 % sample					
3.	120 % sample					
System Stability						
Tailing factor						Not more than 1.8
% Standard RSD					RSD NMT 2.0 %.	
6.0	Robustness:					
6.1	Test	Parameters	Result			
	For Wavelength of UV-Visible Detector:	At 358 nm				
		At 362 nm				
		System Stability				
		Tailing factor				Not more than 1.8
% Standard RSD				RSD NMT 2.0 %.		
6.2	For Flow rate:	At 0.72 ml / min				
		At 0.88 ml / min				
		System Stability				
		Tailing factor			Not more than 1.8	
		% Standard RSD			RSD NMT 2.0 %.	
6.3	For Column temperature:	At 44.5°C				
		At 45.5°C				
		System Stability				
		Tailing factor			Not more than 1.8	
		% Standard RSD			RSD NMT 2.0 %.	
7.0	Solution stability:					



PHARMA DEVILS

QUALITY ASSURANCE DEPARTMENT

ANALYTICAL METHOD VALIDATION REPROT FOR LIVOFLOXACIN TABLETS (FOR ASSAY)

S. No	Analytical parameter		Result			Acceptance criteria
	Parameters	Results (%)	Tailing factor (NMT-1.8)	% Standard RSD (NMT-2.0 %)	Fresh and verification standard difference (NMT-2.0 %)	Recovery value (NMT-2.0 %) from the initial result
	Initial					
	4 hours.					
	8 hours.					
	12 hours.					
	24 hours.					

Conclusion:

Prepared by

Date:

Checked by

Date: