



PHARMA DEVILS
QUALITY ASSURANCE DEPARTMENT

ANNEXURE VII

Batch No. :					
Name of equipment	Results Microbiological contamination by Swabs Technique				
	Location name	Equipment ID No.	Sample ID. No.	Acceptance criteria (100 cm ²)	
				TAMC (NMT 30 cfu)	Pathogen (Should be absent)
SS Scoop	Inner Surface	*	M1	*	*
Sampling Thief	Outer cover inner surface	*	M1	*	*
	Outer surface of inner rod		M2	*	*
Vibrosifter	Inner Surface of Sifter Bowl	*	M1	*	*
	Sifter Sieve surface		M2	*	*
Conta Bin (250 Ltr.)	Inner corner	*	M1	*	*
	Discharge valve of Bin		M2	*	*
Compression M/C	Turret		M1	*	*
	Inner surface of Hopper	*	M2	*	*
	Inside portion of Force feeder		M3	*	*
	Rubber gasket of Y chute side-I		M4	*	*
Tablets Deduster	Inner surface of spiral path of De-duster outlet (charging side)	*	M1	*	*
Metal Detector	Inlet of Metal Detector	*	M1	*	*
Blisters Packing M/C	Hopper	*	M1	*	*
	Rotatory Vibrator		M2	*	*
Deblistering M/c	Inner surface of tablets collecting tray	*	M1	*	*

Checked By (Sign & Date)

Verified By (Sign & Date)