

QUALITY ASSURANCE DEPARTMENT

# TRACEABILITY MATRIX FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

# **TRACEABILITY MATRIX**

# FOR

# **COMPUTER SYSTEM**

# OF

QC (Viscometer)

| System Name    | QC (Viscometer) |
|----------------|-----------------|
| System ID      |                 |
| Location       | Instrument Lab  |
| Effective Date |                 |

Document No.:



Pharma Devils

### QUALITY ASSURANCE DEPARTMENT

# TRACEABILITY MATRIX FOR COMPUTER SYSTEM OF VISCOMETER

System Name: QC (Viscometer)

System ID:

## TABLE OF CONTENTS

| 1.  | PREPARATION AND APPROVALS   | .3  |
|-----|-----------------------------|-----|
| 2.  | REVISION HISTORY            | .4  |
| 3.  | OBJECTIVE AND SCOPE         | .4  |
| 4.  | INTRODUCTION                | .4  |
| 5.  | TRACEABILITY COLUMN DETAILS | .4  |
| 6.  | REFERENCES                  | .5  |
| 7.  | ROLE AND RESPONSIBILITY     | .6  |
| 8.  | TRACEABILITY MATRIX         | .7  |
| 9.  | ABBREVIATION                | . 8 |
| 10. | SUMMARY & CONCLUSION        | . 8 |
|     |                             |     |



QUALITY ASSURANCE DEPARTMENT

### TRACEABILITY MATRIX FOR COMPUTER SYSTEM OF VISCOMETER

# System Name: QC (Viscometer)

System ID:

## 1. PREPARATION AND APPROVALS:

The signature listed below indicates the preparation and approval of this Traceability Matrix. This approval is joint responsibility of listed functional areas.

|             | PROTOCOL DEVELOPMENT | SIGN / DATE |
|-------------|----------------------|-------------|
| Name        | :                    |             |
| Designation | :                    |             |
|             |                      |             |

|             | PROTOCOL REVIEW AND APPROVAL (M/S) |
|-------------|------------------------------------|
| Sign / Date | :                                  |
| Name        | :                                  |
| Designation | :                                  |
|             | Engineering                        |
| Sign / Date | :                                  |
| Name        | :                                  |
| Designation | :                                  |
|             | IT                                 |
| Sign / Date | :                                  |
| Name        | :                                  |
| Designation | :                                  |
|             | Quality Control                    |

|             | PROTOCOL APPROVAL (M/S) |
|-------------|-------------------------|
| Sign / Date | :                       |
| Name        | :                       |
| Designation | :                       |
|             | Quality Assurance       |



# PHARMA DEVILS GUALITY ASSURANCE DEPARTMENT

# TRACEABILITY MATRIX FOR COMPUTER SYSTEM OF VISCOMETER

#### System Name: QC (Viscometer)

System ID:

### 2. REVISION HISTORY:

| Date | Supersedes | Reason for Revision |
|------|------------|---------------------|
|      |            |                     |
|      |            |                     |

### **3. OBJECTIVE AND SCOPE:**

The traceability matrix is prepared for the Computer system Based of QC\_VISCO. Project planned at Metrocraft, Baddi . The traceability matrix contains all the traceability mentioned in system requirement specifications.

This document is dynamic in nature and can be changed during the course of the project with suitable revision number.

### 4. INTRODUCTION:

To provide the relation between system requirement specifications and installation & operational qualification documents to ensure that all the requirement as per the system requirement specifications.

## 5. TRACEABILITY COLUMN DETAILS:

SRS Test Reference No.

System Requirement Specifications requirement should be listed.

SRS Test Description

System Requirement Specifications description should be listed.

 $\blacktriangleright GxP Impact? (Y/N)$ 

If Y, then there must be a test reference in column 5, or a reference showing that this requirement is verified in some other way.

 $\blacktriangleright \quad \text{Other Impact (Y/N)}$ 

The system may require some formal verification or testing for reasons other than GxP, and for which it would be good practice to trace from requirement through design to testing. It is recommended this column notes the reason for other impact.

Verification Test Reference No.



QUALITY ASSURANCE DEPARTMENT

## TRACEABILITY MATRIX FOR COMPUTER SYSTEM OF VISCOMETER

#### System Name: QC (Viscometer)

System ID:

A reference must be present where GxP impact Y.

#### Remarks

Mentioned any remark that add information particularly where reference needs to be made to additional testing or requirements that have arisen as part of the excerise.

### 6. **REFERENCES**:

The publication listed below form part of this protocol's reference documents. Each publication shall be the latest revision in effect on the date this protocol is approved for execution unless noted otherwise. Except as modified by the requirements specified herein or the details of the drawings, work included in this protocol shall conform to the applicable provisions of these publications.

|                                | Good Automated Manufacturing Practices, Version 5,    |  |
|--------------------------------|---|--|
| GAMP 5                         | Guideline document for Automated Systems from         |  |
|                                | International Society of Pharmaceutical Engineering.  |  |
| 21 Code of Federal             | Current Good Manufacturing Practice in Manufacturing, |  |
| <b>Regulations (CFR), Part</b> | Processing, Packing, or Holding off Drugs; General    |  |
| 210                            |   |  |
| 21 Code of Federal             | Current Good Manufacturing Practice for finished      |  |
| <b>Regulations (CFR), Part</b> | Pharmaceuticals                                       |  |
| 211                            |   |  |
| WHO                            | Appendix 5, validation of computerized systems.       |  |
| VP                             | -   |  |
| SRS                            | -   |  |
| RA                             | -   |  |
| IQ                             | -   |  |
| OQ                             | -   |  |
| PQ                             | -   |  |



QUALITY ASSURANCE DEPARTMENT

## TRACEABILITY MATRIX FOR COMPUTER SYSTEM OF VISCOMETER

#### System Name: QC (Viscometer)

System ID:

# 7. ROLE AND RESPONSIBILITY:

The validation team comprising of representative from each of the following departments should be responsible for overall compliance with this validation plan.

| Department        | Responsibilities   |  |  |
|-------------------|--|--|--|
| Validation Agency | $\succ$ To collect data necessary for the generation, execution of this protocol |  |  |
| ()                | from M/s   |  |  |
|                   | > To prepare the traceability matrix.  |  |  |
|                   | $\succ$ To identify the system inventory and document the traceability matrix.   |  |  |
|                   | > Comply with regulatory / Guidelines / Standards / validation plan              |  |  |
|                   | requirements throughout the validation life cycle.                               |  |  |
|                   | > To submit validation documents for approval.                                   |  |  |
| Engineering       | Project Management and Planning.   |  |  |
| ( <b>M</b> /s)    | > To provide the necessary data for qualification activities.                    |  |  |
|                   | To co-ordinate during execution of qualification activities.                     |  |  |
|                   | > To review the validation documents.  |  |  |
| IT                | To provide the necessary data for qualification activities.                      |  |  |
| ( <b>M</b> /s)    | > To co-ordinate during execution of qualification activities.                   |  |  |
|                   | > To review the validation documents.  |  |  |
| Quality Control   | > To provide the necessary data for qualification activities.                    |  |  |
| (M/s)             | > To co-ordinate during execution of qualification activities.                   |  |  |
|                   | > To review the validation documents.  |  |  |
| Quality Assurance |  |  |  |
| (M/s)             | > To approve and authorize the validation documents.                             |  |  |



QUALITY ASSURANCE DEPARTMENT

# TRACEABILITY MATRIX FOR COMPUTER SYSTEM OF VISCOMETER

### System Name: QC (Viscometer)

System ID:

# 8. TRACEABILITY MATRIX:

| SRS<br>Reference<br>No. | SRS Test Description                                       | GxP<br>Impact?<br>(Y/N) | Other<br>Impact?<br>(Y/N) | Verification of<br>Test Reference<br>(IQ & OQ) | Remarks  |
|-------------------------|--|-------------------------|---------------------------|--|----------|
| 10.1                    | Hardware Components  | Y                       | Ν                         |  | Found Ok |
| 10.2                    | Software Components  | Y                       | Ν                         |  | Found Ok |
| 10.3                    | Capacity Requirement                                       | Y                       | N                         |  | Found Ok |
| 10.4                    | Power Utility  | Y                       | N                         |  | Found Ok |
| 10.5                    | Environmental Condition                                    | Y                       | N                         |  | Found Ok |
| 10.6                    | Communication Link<br>Between Server To<br>Computer System | Y                       | N                         |  | Found Ok |
| 10.7                    | Window Security  | Y                       | Ν                         |  | Found Ok |
| 10.8                    | Password Security  | Y                       | Ν                         |  | Found Ok |
| 10.9                    | Verification of User Level<br>and Rights                   | Y                       | N                         |  | Found Ok |
| 10.10                   | System Response Failure                                    | Y                       | Ν                         |  | Found Ok |
| 10.11                   | Electronic Data Security                                   | Y                       | N                         |  | Found Ok |
| 10.12                   | Audit Trail  | Y                       | Ν                         |  | Found Ok |
| 10.13                   | Report Generation  | Y                       | N                         |  | Found Ok |
| 10.14                   | Data Back Up   | Y                       | N                         |  | Found Ok |
| 10.15                   | User Prevented From<br>Alternating Date and Time           | Y                       | N                         |  | Found Ok |
| 10.16                   | 21 CFR part 11 Clauses                                     | Y                       | Ν                         |  | Found Ok |
| 10.17                   | Control Loop Test  | Y                       | N                         |  | Found Ok |

QUALITY ASSURANCE DEPARTMENT

# TRACEABILITY MATRIX FOR COMPUTER SYSTEM OF VISCOMETER

# System Name: QC (Viscometer)

System ID:

### 9. ABBREVIATION:

| Abbreviations | Description                          |  |
|---------------|--------------------------------------|--|
| GMP           | Good Manufacturing Practices         |  |
| IQ            | Installation Qualification           |  |
| OQ            | Operation Qualification              |  |
| PQ            | Performance Qualification            |  |
| QA            | Quality Assurance                    |  |
| ТМ            | Traceability Matrix                  |  |
| SOP           | Standard Operating Procedure         |  |
| SRS           | System Requirement and Specification |  |
| QA            | Quality Assurance                    |  |
| QC            | Quality Control                      |  |
| WHO           | World Health Organization            |  |

### **10. SUMMARY & CONCLUSION:**

Compiled by: \_\_\_\_\_

Date:\_\_\_\_\_

Document No.:



